



Anestand®

Features and Benefits:

- Anestand was designed by anesthesiologists who questioned the practice and safety of using the arm board, shirt pocket, or patient's chest to organize and hold critical supplies.
- Anestand can be positioned and repositioned in multiple configurations creating an adaptable workspace in multiple anesthetic situations.
- The customized tray holds items for induction or procedures.
- The IV clips can be used to keep your IV port immediately available and clean and designed to assist in complying with the JCAHO "Scrub the Hub!" initiative.
- Corner grips can be used to organize ETT and suction tubing.
- Circuit clips assist in supporting airway circuit.

Assembly/Operation:

- Anestand is assembled from 2 main parts: the tray and flexible support post with Gibbs clamp.
- Anestand is assembled by sliding the tray clip into the tray clip receptacle on the underside of the tray. (Figure 1)
- Anestand can easily be mounted to an IV pole or OR bed rail using the versatile Gibbs clamp. (Figures 2 and 3)




Figure 1




Figure 2




Figure 3



Warnings:

- Anestand is NOT MRI compatible.
- Do NOT autoclave. Heat sterilization will result in damage to Anestand as well as sterilization equipment.
- Anestand was not designed for direct contact with the patient. Position Anestand in a fashion to avoid direct patient contact.
- Care is needed to avoid dropping equipment when Anestand is repositioned.
- Anestand was not designed to hold more than 2-lbs when the flexible support post is straight or more than 1-lb when the support post is positioned at a 90 degree angle.
- If any part of Anestand is damaged, or if the tray cannot be fully cleaned, it must be replaced. Replacement trays are available. Please visit anestand.com
- As with any circuit stand, it is important to keep your Anestand anesthesia stand clear of the patient. This is especially true of an emerging patient. Anestand recommends that your anesthesia stand not be fixed to the bed on emergence, but rather somewhere clear of the patient (like the IV pole for example).

Cleaning and Maintenance:

- All standard disinfectants, disinfectant wipes, mild soap and water solutions, and mild detergents are safe to use in cleaning Anestand.
- Use a sterilization protocol that complies with your facility's infection control and risk management policies for anesthesia machine surfaces.
- Anestand should be inspected before every use to ensure that the Gibbs clamp securely closes, the flexible support post is functional, and all parts are securely fashioned.
- All surfaces must be cleaned after each use.
- IV clips, corner grips, and the tray clip receptacle must be free of any debris and be cleaned with facility approved cleaning materials

Technical Specification:

- Latex Free
- Tray made of 100% non porous polypropylene
- Made in the USA
- Support post covered in medical grade non porous silicone
- FDA Registered
- US Patent Pending: 15/789,979
- European patent office application number: 17862480.5
- Canadian patent application number : 3,038,004
- Please visit anestand.com for ordering information