

# Patient Drainage Bags



**M04200** Drainage Bag  
2000mL



**M08101** Drainage Bag  
4000mL

- **120cm Detachable Drain Line with Rotating Male Luer Fitting**  
Direct and simple hook up to drainage catheters offered by MedQuest or any standard female luer connection
- **Roller Clamp Assembly**  
Located on the drain line for quick and easy flow rate control for patients experiencing discomfort during drainage
- **Built-in Rotating Hanger, Carrying Handle, and Tie Cord**  
Increases ease of use and handling in a variety of settings
- **Integrated Drip Chamber and Anti-Reflux Valve**  
Prevents fluid from returning to the patient and decreases risk of potential ascending line infections
- **Bacteria-Proof and Liquid-Tight Ventilation**  
Ensures reliable flow performance and optimal pressure balance inside the bag
- **Outlet Housing to Secure Outlet Drainage Tubing**  
Guards against contamination of surroundings and increases floor clearance
- **Non-Drip Bottom Outlet with Lever Clamp**  
Enables hygienic and rapid emptying of bag

© 2016 by MedQuest Medical, Inc. All Rights Reserved.