

# Medquest Universal Drainage Kit S55148KT (1000ml)

The Medquest Universal Drainage Kit gives Hospitals and Home Care Providers the ability to drain all commonly available tunneled Catheters for Pleural Effusions and Malignant Ascites.



- The SA5200 Drain Line Adapter may be also used to access and drain all non-leur catheter valves requiring a stylet access tip.

- We have included a Universal Catheter Cap in the Kit (PCA-1050) for catheters requiring protection between dressing changes.



- The S55148, 1000ml bottle, can be used to drain all catheters with female leur lock fittings.

**Medquest Complete Dressing Kit**  
(Product Code 2976C)  
**may be ordered separately.**

- Kit Contents Include: 1 each S55148 (1000ml) Sterile Bottle with attached Drain Line.  
1 each SA5200 Safety Drain Line Adapter  
1 each PCA-1050 Universal Catheter Cap