



Pleural & Peritoneal Drainage Systems

To place an order or to learn more about how MedQuest Medical can assist in providing additional information on this valuable treatment option, contact:

MedQuest Customer Service

Email: customer.service@medquestmedical.com

Phone: 888-864-4442

Rev 02/2018



Distributed exclusively by:
MedQuest Medical Inc.,
www.medquestmedical.com | Ph./Fax : 888 864-4442



Home management of recurrent pleural effusions & malignant ascites can now be accomplished by using MedQuest's cost effective vacuum bottles or patient drainage bags.

Why choose a long-term indwelling catheter featuring a standard luer connection?

- Secure connection
- Nurse Efficiency
- Cost-Effective Supplies
- Replaceable Valve
- Ease of use

MedQuest Medical Inc

**Pleural / Peritoneal
CATHETER DRAINAGE
SYSTEMS**

Why choose MedQuest drainage options?

Versatile

- Our product offering consists of vacuum bottle & drainage bag options that allow the physician to match your patient's clinical needs with our integrated collection solutions

User friendly, well-tolerated

- Drainage can be easily managed through homecare nurse providers, compliant patients & family members.
- Physicians, nurses and patients can be confident that this therapy is successfully used by leading hospitals throughout the country.

Invoke positive change

- You can now choose the cost-effective option of a direct connection to 2000mL or 4000mL patient drainage bags for outpatient drainage of malignant ascites.
- At home drainage reduces hospital length of stay(LOS), thereby minimizing the risk of healthcare-associated infections (HAI) for the immunocompromised patient.
- Reduced LOS & risk of HAI reduce costs to the overall health care system.

Education

- Extensive, ongoing nurse and patient education is a priority to MedQuest Medical.
- Comprehensive nurse and patient educational material is available upon request.

Product Ordering Information

<i>Product No.</i>	<i>Description</i>	<i>Quantity</i>
S55120	600mL Drainage Bottle, Sterile	10/box
S55148	1000mL Drainage Bottle, Sterile	10/box
S55149	1000mL Non-Sterile Drainage Bottle	10/box
SA5200	Safety Drain Line Adapter, Sterile	50/box
DCS-100	Pleural Catheter Kit, Sterile	Each
DCS-101	Peritoneal Catheter Kit, Sterile	Each
2976	Dressing Kit, Sterile	10/box
797	Luer Cap, Sterile	25/box
DRCV-201	Drainage Catheter Valve, Sterile	5/box
M04200	2000mL Drainage w/Drip Chamber & Drain Line	Each
M08101	4000mL Drainage w/Drip Chamber & Drain Line	Each
S55300	Universal Valve Cap	Each

Patient Drainage Bags



M04200 Drainage Bag
2000mL



M08101 Drainage Bag
4000mL

- **120cm Detachable Drain Line with Rotating Male Luer Fitting**
Direct and simple hook up to drainage catheters offered by MedQuest or any standard female luer connection
- **Roller Clamp Assembly**
Located on the drain line for quick and easy flow rate control for patients experiencing discomfort during drainage
- **Built-in Rotating Hanger, Carrying Handle, and Tie Cord**
Increases ease of use and handling in a variety of settings
- **Integrated Drip Chamber and Anti-Reflux Valve**
Prevents fluid from returning to the patient and decreases risk of potential ascending line infections
- **Bacteria-Proof and Liquid-Tight Ventilation**
Ensures reliable flow performance and optimal pressure balance inside the bag
- **Outlet Housing to Secure Outlet Drainage Tubing**
Guards against contamination of surroundings and increases floor clearance
- **Non-Drip Bottom Outlet with Lever Clamp**
Enables hygienic and rapid emptying of bag

System Features and Benefits

1) Vacuum Bottle

PreVac 600mL and 1000mL bottles effectively drain effusions when required using active negative pressure.

2) Vacuum Indicator

Green vacuum indicator confirms negative pressure is available for drainage.

3) Drain Line with Roller Clamp

120cm long drain line provides options for bottle placement around procedure area while draining. Roller clamp puts ease of control in hands of the caregiver for increased patient comfort.

4) Standard Luer Safety Valve

Standard female luer connection on safety valve permits secure attachment to the male luer end of tubing set for consistent and convenient drainage procedures.

5) Polyester Cuff

Cuff promotes subcutaneous tissue ingrowth, helping to reduce risk of infection and securely hold catheter in place.

6) 15.5 Fr Silicone Catheter

Large fenestrations along the soft and flexible catheter promote effective fluid drainage.

7) Safety Drain Line Adapter

Enables connection of MedQuest vacuum bottle to competitive catheter systems.

8) Luer Cap

Application of cap to drainage bottle creates a closed system for safe disposal in outpatient setting. Catheter luer cap to prevent integrity of valve & reduce risk of infection.

9) Drainage Catheter Valve

Ability to replace valve eliminates need to subject patient to additional procedure. Reduced hospital re-admissions & costs associated with catheter related issues.



Comprehensive Catheter Kits

MedQuest Medical Inc. offers catheter kits for both pleural (DCS-100) and peritoneal (DCS-101) applications. Each catheter kit contains all necessary components for catheter insertion.

- 1 J-Tip Guidewire, Graduated
- 1 10cc Syringe L.L.
- 1 20cc Syringe L.L.
- 6 Gauze Pad 4"x4"
- 1 Fenestrated Drape 60"x 43" w/ 10"x6" fenestration
- 1 5-in-1 Drainage Line Adapter
- 1 Drainage Line
- 1 Catheter Stylet
- 1 Drainage Catheter

- 1 18GA Introducer Needle
- 1 16FR Peel-Away Introducer
- 1 Tunneler
- 1 Safety Scalpel #11
- 1 12FR Dilator x 8"
- 1 18GA x 1.5" Safety Needle
- 1 22GA x 1.5" Safety Needle
- 1 4"x4" Split Foam Dressing
- 1 6"x8" Self-Adhesive Dressing Tegaderm



Dressing Kit

MedQuest Medical offers patient dressing kits(2976C) that are sterile packed and complete with all necessary components required to perform dressing changes on patients with indwelling catheters in hospital, home or palliative care setting.

Each dressing kit, wrapped in a CSR wrap, contains:

- 1 ea. 3M Cavilon Barrier Film
- 1 ea. 3M Tegaderm Dressing
- 1 ea. Foam Pad with 2" Split, 4" x 4" square, 3/16" thick
- 8 ea. Alcohol Wipes
- 1 ea. Emergency Tube Pinch Clamp
- 12 ea. Gauze Pads, 4" x 4", 8-ply
- 1 ea. Luer Cap, light blue (inside a small plastic bag)
- 2 pairs Chemo Rated Nitrile Gloves, Medium

