



As soon as  
you use it,  
you'll wonder how  
you ever managed  
without it.



## Presenting Anestand.



*Anestand* is the first and only anesthesia stand designed to:


- Safely, securely and conveniently keep all anesthesia supplies where they're needed.
- Deliver the highest levels of safety, cleanliness, utility and organization.
- End instruments rolling off a patients chest and IV tubing being caught in the wheels of the IV pole.



The many benefits  
of Anestand are  
obvious, every time  
you use it.

**Extremely easy to place and set up.**

The patent pending clamp quickly and easily attaches *Anestand* to the IV pole or OR bed rail with just one hand — no knobs or screws to tighten.




**The recessed tray safely holds and organizes supplies for most procedures and can be used for induction of anesthesia.**

Tray holds: central line kit, nerve block supplies, and spinal/epidural kit. Organized workspace to make an efficient induction.

**Clips keep IV tubing securely in place, so it's easy to "scrub the hub."**

*Anestand* makes it easy to comply with JCAHO's and CDC's "scrub the hub" program. And tubing can't fall onto the OR floor.



**Corners hold suction tubing or suction tip.**

Easily organize EKG, pulse oximetry and blood pressure cords. And all cords stay off the floor.

**The anesthesia circuit clips keep your circuit fixed and firmly in place.**

When the table is in Trendelenburg or reverse Trendelenburg, your circuit remains fixed in place.

**The flexible support post lets you position it just where you need it.**

You can rotate or bend the post, making *Anestand* adaptable to any anesthetic situation.

**It's like an extra pair of hands, wherever and whenever you need them.**

*Anestand* frees up anesthesia technicians for other valuable duties — and works as helping hands.

ANESTAND INDUCTION TRAY #77001  
ANESTAND CLAMP STAND #77002

 **Anestand**<sup>®</sup>  
The anesthesia stand that  
keeps everything right at hand.

(800) 768-7061 info@anestand.com www.anestand.com