



Memory Staple

A body heat activated staple for continuous compression.

The Memory Staple is a Nitinol memory-alloy staple which is activated at patient body temperature. It features a unique S-bend bridge that allowing for even compression across the fusion site while maintaining a low profile against the bone.



Open State



Closed State



Memory Staple

S-Bend Design

The S-Bend bridge of the Memory Staple ensures even compression across the fusion site, while maintaining a low profile against the bone.




Continuous Compression

Manufactured from Nitinol and activated at patient body temperature, the Memory Staple provides dynamic compression during the healing process.



More Sizes = More Flexibility

The Memory Staple is one of the most comprehensive compression staple systems on the market with 17 standard sizes.



Bridge Width	Leg Length	Wire Size
7mm	5mm	1.2mm
7mm	5mm	1.5mm
9mm	7mm	1.2mm
9mm	7mm	1.5mm
11mm	8mm	1.5mm
11mm	10mm	1.5mm
11mm	15mm/13mm	1.5mm
11mm	17mm/15mm	1.5mm
11mm	19mm/17mm	1.5mm
13mm	10mm	1.5mm
15mm	12mm	1.5mm
15mm	12mm	2.0mm
18mm	18mm/15mm	2.0mm
20mm	20mm	2.0mm
20mm	20mm	2.0x3.0mm
25mm	22mm	2.0x3.0mm
30mm	30mm	2.0x3.0mm


Suggested Procedures

Foot and Ankle:

- Hallux IP Fusion: 9x7, 11x8, 11x10, 13x10
- Hallux MPJ Fusion: 13x10, 15x12, 18x18x15
- Lesser PIPJ Fusion: 7x5, 9x7
- Austin, Youngswick, Reverdin Bunionectomies: 11x10, 13x10, 11x15x13
- Base Wedge Osteotomies: 11x10, 13x10, 15x12, 18x18x15
- Akin Osteotomies: 7x5, 9x7
- Cotton Procedure: 20x20, 18x18x15
- Lapidus Procedure (Dorsal): 20x20, 18x18x15
- Lapidus Procedure (Medial): 15x12
- Talonavicular Fusion: 18x18x15, 20x20
- Calcaneal Cuboid Fusion: 15x12, 18x18x15, 20x20
- Dwyer Osteotomy: 20x20

Hand and Wrist:

- DIP Fusion: 7x5, 9x7
- PIP Fusion: 7x5, 9x7
- MCP Fusion: 11x8, 11x10, 13x10
- Capitate/Scaphoid Fusion: 11x8, 11x10, 13x10
- Capitate/Lunate Fusion: 11x8, 11x10, 13x10
- CMC Thumb Fusion: 11x8, 11x10, 13x10
- Third Metacarpal/Capitate Fusion: 11x8, 11x10, 13x10
- Fifth Metacarpal/Hamate Fusion: 11x8, 11x10, 13x10



Visit www.bioproimplants.com or call us at (810) 982-7777



+1-810-982-7777
 info@bioproimplants.com
 2929 Lapeer Road
 Port Huron, MI 48060

