


Technical Specification

Chambers	3
Compression mode	4
Time range	1min-24hours or continuous
Pressure range	40-150mmHg
Hold time	0 - 10s (Step by 1s)
Interval time	10 - 60s (Step by 1s)
Interface	5" touch screen
Rechargeable battery	Yes (Lithium battery : 11.1V 9000mAh)
Adapter	Yes
Input voltage	AC 100 - 240 V 50 - 60 Hz
Dimension	250 (L) x 184 (W) x 150 (H) mm
Weight	2.5 kg

Certification



Longest

Guangzhou Longest Science & Technology Co., Ltd.
 Address: 5&6F, Building B4, No.11, Kaiyuan Avenue, Science City,
 Guangzhou Hi-tech Industrial Development Zone, 510530 Guangzhou, China
 Tel: 86-20-66353999 Fax: 86-20-66353920
 E-mail: export@longest.cn Website: www.longest.cn www.gzlongest.com



LGT-2200DVT

Compression Therapy Device



Health Life Longest Care


Longest



LGT-2200DVT is a special Pressotherapy system that is designed to decrease the risk of DVT(Deep Venous Thrombosis) and Pulmonary Embolism(PE)before, during and after surgery

4 compression modes with reusable & disposable sleeves meet all your need

⊕ Compression Modes



⊕ For Disposable Sleeve



*Notice: For one - leg treatment, please use single - head hose.


⊕ For Reusable Sleeve



*Notice: For two - leg treatment, please use double - head hose.

⊕ Introduction of deep vein thrombosis (DVT)

⊕ Reason



Deep vein thrombosis is one of the most common complications of bed rest, prevention is urgently needed

⊕ Features

- ⊕ 4 Compression Modes
- ⊕ 5" Touch Screen
- ⊕ Reusable & Disposable Sleeves
- ⊕ Gradient Pressure Setting
- ⊕ Intermittent Pneumatic Compression
- ⊕ Pressure of Every Chamber Adjustable
- ⊕ Rechargeable Battery
- ⊕ Powered by AC or DC

